





G06		Vise grips	
		www.grip.de/pdf/THS341_Datasheet.pdf	
		THS341	
	THS341	Small vise grip with chain, 100 N Clamping surface with integrated pyramid jaws 5x6 mm (h x w) Price includes: 1 piece	
	THS341-F	Like THS341, but without chain The coupling is direct mounted on the grip Price includes: 1 piece	
	THS341-Ni THS341-F-Ni	Like THS341 and THS341-xx-F-Ni, but nickel-plated It is suitable for a temperature range from -70°C to 280°C Price includes: 1 piece	
	THS341-10-Ko	Small vise grip with pistons, 100 N Clamping surface with integrated pyramid jaws 5x10 mm (h x w) Price includes: 1 piece	
  		www.grip.de/pdf/TH140_Datasheet.pdf	
		TH140k	
	TH140k	Small vise grip incl. pyramid jaws, 100 N Clamping surface: 15x15 mm Opening: 0-8 mm Weight per piece: 63 gr. Price includes: 1 pair incl. pyramid jaws	
	TH140k+Ko	Pneumatic Version of TH140k like TH140k including pyramid jaws with 1rod Clamping force: 140 N Price includes: 1 pair with pyramid jaws	
	TH140k+2Ko	Pneumatic Version of TH140k like TH140k including pyramid jaws with 2 rods Clamping force: 140 N Price includes: 1 pair with pyramid jaws	
	TH140s	Vertical version of TH140k with Adapter 8 and 15.9 incl. pyramid jaws Clamping surface: 15x15 mm Opening: 0-8 mm Weight per piece: ~45 gr. Price includes: 1 pair with pyramid jaws	
			Additional jaw faces for TH140 Price includes: 1 set = 4 jaw faces
	TH140k-B	Blank (smooth) jaws without coating, 15x15 mm (h x w), 1 set	
	TH140k-BG	Rubber jaws, 15x15 mm (h x w), 1 set	
	TH140k-BP	Pyramid jaws, 15x15 mm (h x w), 1 set (standard jaws)	
 		www.grip.de/pdf/TH227_Datasheet.pdf	
		TH227	
	TH227	Small vise grip, 100 N Aluminium, 60 gram, M5 coupling Opening: till 8 mm (depending on jaws) Price includes: 1 pair (without jaws)	
	TH227-S46	Small vise grip, 100 N Aluminium, 98 gram, M5 coupling Opening: till 46 mm (depending on jaws) Price includes: 1 pair (without jaws)	
			Jaws for TH227 1 Set = 4 jaws Further sizes on request Price includes: 1 Set = 4 jaws
	TH227-B25	Blank (smooth) jaws without coating, 10x25 mm (h x w), 1 set	
TH227-BG25	Rubber jaws 10x25 mm (h x w), 1 set		

	TH227-BP25	Pyramid jaws (0.8x45), 10x25 mm (h x w), 1 set
	TH227-BD25	Diamond coated jaws 10x25 mm (h x w), 1 set
	TH227-B50	Blank (smooth) jaws without coating, 10x50 mm (h x w), 1 set
	TH227-BG50	Rubber jaws 10x50 mm (h x w), 1 set
	TH227-BP50	Pyramid jaws (0.8x45) 10x50 mm (h x w), 1 set
	TH227-BD50	Diamond coated jaws 10x50 mm (h x w), 1 set
	TH227-B100	Blank (smooth) jaws without coating, 10x100 mm (h x w), 1 set
	TH227-BG100	Rubber jaws 10x100 mm (h x w), 1 set
	TH227-BP100	Pyramid jaws (0.8x45) 10x100 mm (h x w), 1 set
	TH227-BD100	Diamond coated jaws 10x100 mm (h x w), 1 set
THS637-159-Km5	Female Adapter (see THS637-159-Km5) Price for 1 pair (2 adapters)	
TH9		
	TH9	Small screw action grip for thin wires or ropes, 2.5 kN Opening: 0-5 mm Body: aluminium Weight per piece: ~ 380 gr Price includes: 1 Pair of grips and carriers for inset jaws (but without inset jaws)
	TH9+Ko	Pneumatic version of TH9 for wires, 2.5 kN Opening: 0-5 mm Clamping force: 1.3 kN Body: aluminium Weight per piece: ~ 550 gr Price includes: 1 Pair of grips and carriers for inset jaws (but without inset jaws)
	TH9+2Ko	Pneumatic version of TH9 with 2 pistons for wires, 2.5 kN Opening: 0-5 mm Clamping force: 1.3 kN Body: aluminium Weight per piece: ~ 710 gr Price includes: 1 Pair of grips and carriers for inset jaws (but without inset jaws)
	Inset jaw faces for TH9 - 1 set = 4 jaw faces	
	TH9-EB	Blank (smooth) inset jaws 37x15 mm (h x w), 1 set
	TH9-EBD	Diamond coated inset jaws 37x15 mm (h x w), 1 set
http://www.grip.de/g06/TH240k/		
www.grip.de/pdf/TH240k_Datasheet.pdf		
TH240k		
 	TH240k	Vise grip U-Form, 2.5 kN Opening: 0-10 mm Price includes: 1 pair (without jaws)
	TH240k-S20	Vise grip U-Form, 2 kN Opening: 0-20 mm Price includes: 1 pair (without jaws)
	TH240k-S30	Vise grip U-Form, 1 kN Opening: 0-30 mm Price includes: 1 pair (without jaws)
	TH240k-S50	Vise grip U-Form, 1 kN Opening: 0-50 mm Price includes: 1 pair (without jaws)
Pneumatic Versions of TH240k		

TH240k+Ko	Vise action grip U-Form with 1 pneumatic rod Ø 50 mm, 1.2 kN Opening: till 10 mm (depending on jaw faces) Price includes: 1 pair (without jaws)
TH240k+2Ko	Vise action grip U-Form with 2 pneumatic rods Ø 50 mm, 1.2 kN Opening: till 10 mm (depending on jaw faces) Price includes 1 pair (without jaws)
TH240k-S20+Ko	Vise action grip U-Form with 1 pneumatic rod Ø 50 mm, 1.2 kN Opening: till 20 mm Price includes: 1 pair (without jaws)
TH240k-S20+2Ko	Vise action grip U-Form with 2 pneumatic rods Ø 50 mm, 1.2 kN Opening: till 20 mm Price includes: 1 pair (without jaws)
TH240k-S30+Ko	Vise action grip U-Form with 1 pneumatic rod Ø 50 mm, 1.2 kN Opening: till 30 mm Price includes: 1 pair (without jaws)
TH240k-S30+2Ko	Vise action grip U-Form with 2 pneumatic rods Ø 50 mm, 1.2 kN Opening: till 30 mm Price includes: 1 pair (without jaws)
TH240k-S50+Ko	Vise action grip U-Form with 1 pneumatic rod Ø 50 mm, 1.2 kN Opening: 0 - 50 mm Price includes: 1 pair (without jaws)
TH240k-S50+2Ko	Vise action grip U-Form with 2 pneumatic rods Ø 50 mm, 1.2 kN Opening: 0 - 50 mm Price includes: 1 pair (without jaws)
	Jaw faces for TH240k Price includes: 1 Set = 4 jaws
TH240k-B	Blank (smooth) jaws without coating, 30x30 mm, 1 set
TH240k-BG	Rubber jaws, 30x30 mm, 1 set
TH240k-BP	Pyramid (serrated) jaws, 30x30 mm, 1 set
TH240k-BV2	V-jaws, 30x30 mm, 1 set
TH240k-BV3	V-jaws, 30x30 mm, 1 set
TH240k-BW	Wave jaws, 30x30 mm, 1 set
TH240k-BD	Diamond coated jaws, 30x30 mm, 1 set
TH240k-BLC	Line contact jaws, radius 3, 30x30 mm Price includes: 1 set including 2xBG rubber jaws and 2xBLC line contact jaws
TH240k-B50	Blank (smooth) jaws without coating, 30x50 mm (h x w), 1 set
TH240k-BG50	Rubber jaws, 30x50 mm (h x w), 1 set
TH240k-BP50	Pyramid (serrated) jaws, 30x50 mm (h x w), 1 set
TH240k-BW50	Wave jaws, 30x50 mm (h x w), 1 set
TH240k-BD50	Diamond coated jaws, 30x50 mm (h x w), 1 set
TH240k-B80	Blank (smooth) jaws without coating, 30x80 mm (h x w), 1 set
TH240k-BG80	Rubber jaws, 30x80 mm (h x w), 1 set
TH240k-BP80	Pyramid (serrated) jaws, 30x80 mm (h x w), 1 set
TH240k-BW80	Wave jaws, 30x80 mm (h x w), 1 set
TH240k-BD80	Diamond coated jaws, 30x80 mm (h x w), 1 set
TH240k-B100	Blank (smooth) jaws without coating, 30x100 mm (h x w), 1 set
TH240k-BG100	Rubber jaws, 30x100 mm (h x w), 1 set
TH240k-BP100	Pyramid (serrated) jaws, 30x100 mm (h x w), 1 set
TH240k-BW100	Wave jaws, 30x100 mm (h x w), 1 set
TH240k-BD100	Diamond coated jaws, 30x100 mm (h x w), 1 set


see also G12 for pneumatic

		Inset jaw faces for TH240k Quick change of jaw surface Price includes: 1 Set = 4 jaws
	TH240k-GT	Carrier for inset jaws Price includes: 1 Set (4 carriers)
	TH240k-EBB	Blank (smooth) insert jaws blank, 30x36 mm (h x w), 1 set
	TH240k-EBG	Rubber coated inset jaws, 30x36 mm (h x w), 1 set
	TH240k-EBP	Pyramid (serrated) inset jaws, 30x36 mm (h x w), 1 set
	TH240k-EBD	Diamond coated inset jaws, 30x36 mm (h x w), 1 set
	TH240k-EBG50	Rubber coated inset jaws, 50x36 mm (h x w), 1 set
	TH240k-EBP50	Pyramid (serrated) inset jaws, 50x36 mm (h x w), 1 set
	TH240k-EBD50	Diamond coated inset jaws, 50x36 mm (h x w), 1 set
	TH240k-EBD50	Wave jaws, 50x36 mm (h x w), 1 set
		www.grip.de/pdf/TH240g_Datasheet.doc
		TH240g
	TH240g	Vise grip U-Form Tensile force: 5 kN Opening: till 30 mm (depending on jaw faces) Price includes: 1 pair (without jaws)
	TH240g-S50	Vise grip U-Form Tensile force: 5 kN Opening: till 50 mm (depending on jaw faces) Price includes: 1 pair (without jaws)
	TH240g-S100	Vise grip U-Form Tensile force: 5 kN Opening: 10-100 mm (depending on jaw faces) Price includes: 1 pair (without jaws)
		Pneumatic Versions www.grip.de/pdf/TH240g+Ko_Datasheet.pdf
	TH240g+Ko	Vise action grip U-Form with 1 pneumatic rod Ø 70 mm Clamping force: 2.6 kN Opening: till 30 mm (depending on jaw faces) Price includes: 1 pair (without jaws)
	TH240g+2Ko	Vise action grip U-Form with 2 pneumatic rods Ø 70 mm Clamping force: 2.6 kN Opening: till 30 mm (depending on jaw faces) Price includes: 1 pair (without jaws)
	TH240g-el	Quick release lever (incl. remodelling), 1 pair
		Jaw Faces for TH240g + TH240g-S50 Price includes: 1 Set = 4 jaws
	TH240g-B	Blank (smooth) jaws without coating, 30x50 mm (h x w), 1 set
	TH240g-BG	Rubber jaws 30x50 mm (h x w), 1 set
	TH240g-BP	Pyramid jaws 30x50 mm (h x w), 1 set
	TH240g-BV2	V-jaws for round sample Ø 2-15 mm, 30x50 mm (h x w), 1 set
	TH240g-BV5	V-jaws for round sample Ø 5-30 mm, 30x50 mm (h x w), 1 set
	TH240g-BW	Wave jaws 30x50 mm (h x w), 1 set
	TH240g-BD	Diamond coated jaws 30x50 mm (h x w), 1 set
		Jaw Faces for TH240g + TH240g-S50 + TH240g-S100 Price includes: 1 Set = 4 jaws
	TH240g-B80	Blank (smooth) jaws without coating, 30x80 mm (h x w), 1 set
	TH240g-BG80	Rubber jaws 30x80 mm (h x w), 1 set
	TH240g-BP80	Pyramid jaws 30x80 mm (h x w), 1 set

siehe auch G12



TH240g-BW80	Wave jaws 30x80 mm (h x w), 1 set
TH240g-BD80	Diamond coated jaws 30x80 mm (h x w), 1 set
TH240g-B100	Blank (smooth) jaws without coating, 30x100 mm (h x w), 1 set
TH240g-BG100	Rubber jaws 30x100 mm (h x w), 1 set
TH240g-BP100	Pyramid jaws 30x100 mm (h x w), 1 set
TH240g-BW100	Wave jaws 30x100 mm (h x w), 1 set
TH240g-BD100	Diamond coated jaws 30x100 mm (h x w), 1 set
	Inset jaw faces - quick change of jaw surface, for TH240g + TH240g-S50 Price includes: 1 Set = 4 jaws
TH240g-GT	Carriers for inset jaws, 1 set (4 carriers)
TH240g-EBB	Blank (smooth) inset jaws, 30x50 mm (h x w), 1 set
TH240g-EBG	Rubber inset jaws, 30x50 mm (h x w), 1 set
TH240g-EBP	Pyramid inset jaws, 30x50 mm (h x w), 1 set
TH240g-EBD	Diamond coated inset jaws, 30x50 mm (h x w), 1 set



www.grip.de/g06/TH56-k/



	www.grip.de/pdf/TH56k_Datasheet.pdf
	TH56k
TH56k	Vise grip Tensile force: 1 kN Price includes: 1 pair of grips (without jaws)
	Pneumatic version of TH56k (see G12)
TH56k+Ko	Vise action grip with 1 pneumatic rod Ø 50 mm Clamping force: 1.2 kN Price includes: 1 pair (without jaws)
TH56k+2Ko	Vise action grip with 2 pneumatic rods Ø 50 mm Clamping force: 1.2 kN Price includes: 1 pair (without jaws)
TH56k-B..	Jaw faces for TH56k Price includes: 1 set (4 jaws)
TH56k-B	Blank (smooth) jaws without coating, 1 set
TH56k-BG	Rubber jaws, 1 set
TH56k-BP	Pyramid (serrated) jaws, 1 set
TH56k-BV2	V-jaws for round samples Ø 2-10 mm, 1 set
TH56k-BV3	V-jaws for round samples Ø 3-10 mm, 1 set
TH56k-BW	Wave jaw faces for textiles, 1 set
TH56k-BD	Diamond coated jaw faces, 1 set
TH56k-E...	Inset jaws for TH56k Price includes: 1 set (4 inset jaws)
TH56k-GT	Carrier for inset jaws, 1 set
TH56k-EBB	Blank (smooth) inset jaws without coating, 1 set
TH56k-EBG	Rubber coated inset jaws, 1 set
TH56k-EBP	Pyramid inset jaws, 1 set
TH56k-EBD	Diamond coated inset jaws, 1 set

siehe auch G12



	www.grip.de/pdf/TH56g_Datasheet.pdf
	TH56g
TH56g	Vise grip Tensile force: 5 kN Price includes: 1 pair (without jaws)



TH56g-H200	Vise grip Height of the grip: 220 mm Tensile force: 5 kN Price includes: 1 pair (without jaws)
-------------------	---

[Pneumatic version www.grip.de/pdf/TH56g_Datasheet.pdf](http://www.grip.de/pdf/TH56g_Datasheet.pdf)



TH56g+Ko	Pneumatic grip with 1 rod Ø 70 mm Clamping force: 2.6 kN Price includes: 1 pair (without jaws)
-----------------	--

TH56g+2Ko	Pneumatic grip with 2 pneumatic rods Ø 70 mm Clamping force: 2.6 kN Price includes: 1 pair (without jaws)
------------------	---

TH56g-H200+ko	Pneumatic grip with 1 pneumatic rod Ø 70 mm Height of the grip: 220 mm Clamping force: 2.6 kN Price includes: 1 pair (without jaws)
----------------------	--

TH56g-H200+2ko	Pneumatic grip with 1 pneumatic rod Ø 70 mm Height of the grip: 220 mm Clamping force: 2.6 kN Price includes: 1 pair (without jaws)
-----------------------	--



TH56g-el	Quick release lever incl. remodelling Optional for all TH56g grips Surcharge for 2 levers
-----------------	---

Jaw faces for TH56g
Price includes: 1 Set = 4 jaw faces



TH56g-B	Blank (smooth) jaws without coating, 30x50 mm (h x w), 1 set
----------------	--

TH56g-BG	Rubber jaws 30x50 mm (h x w), 1 set
-----------------	-------------------------------------

TH56g-BP	Pyramid (serrated) jaws 30x50 mm (h x w), 1 set
-----------------	---

TH56g-BV2	V-jaws for round samples Ø 2-15 mm, 30x50 mm (h x w), 1 set
------------------	---

TH56g-BV5	V-jaws for round samples Ø 5-11 mm, 30x50 mm (h x w), 1 set
------------------	---

TH56g-BW	Wave jaws 30x50 mm (h x w), 1 set
-----------------	-----------------------------------

TH56g-BD	Diamond coated jaws 30x50 mm (h x w), 1 set
-----------------	---

Inset jaw faces
Quick change of jaw surface
Price includes 1 Set = 4 jaw faces



TH56g-GT	Carrier for inset jaws, 1 set (4 carriers)
-----------------	--

TH56g-EBB	Blank (smooth) inset jaws, 30x50 mm (h x w), 1 set
------------------	--

TH56g-EBG	Rubber inset jaws, 30x50 mm (h x w), 1 set
------------------	--

TH56g-EBP	Pyramid inset jaws, 30x50 mm (h x w), 1 set
------------------	---

TH56g-EBD	Diamond coated inset jaws, 30x50 mm (h x w), 1 set
------------------	--

<http://www.grip.de/q06/TH154/>

www.grip.de/pdf/TH154_Datasheet.pdf

TH154



TH154-10	Vise grip, aluminium Tensile force: 10kN, weight: 2.36 kg Opening: till 34 mm (depending on the type of jaws) Price includes: 1 pair (without jaws)
-----------------	--

TH154-20	Vise grip, steel Tensile force: 20 kN, weight: 5 kg Opening: till 34 mm (depending on the type of jaws) Price includes: 1 pair (without jaws)
-----------------	--

Jaw faces for TH154
Price includes: 1 Set = 4 jaw faces

TH154-B	Blank (smooth) jaws without coating, 40x60 mm (h x w), 1 set
----------------	--

TH154-BG	Rubber jaws 40x60 mm (h x w), 1 set
-----------------	-------------------------------------

TH154-BP	Pyramid (serrated) jaws 40x60 mm (h x w), 1 set
-----------------	---

TH154-BV4	V-jaws for round samples Ø 4-30 mm, 40x60 mm (h x w), 1 set
------------------	---



TH154-BV8	V-jaws for round samples Ø 8-34 mm, 40x60 mm (h x w), 1 set
TH154-BW	Wave jaws 40x60 mm (h x w), 1 set
TH154-BD	Diamond coated jaws 40x60 mm (h x w), 1 set
TH154-BWC	Wave jaws for test of bitumen (with water cooling) 40x60 mm (h x w), 1 set incl. tubes and fittings
TH154-B100	Blank (smooth) jaws without coating, 40x100 mm (h x w), 1 set
TH154-BG100	Rubber jaws, 40x100 mm (h x w), 1 set
TH154-BP100	Pyramid (serrated) jaws, 40x100 mm (h x w), 1 set
TH154-BW100	Wave jaws, 40x100 mm (h x w), 1 set
TH154-BD100	Diamond coated jaws, 40x100 mm (h x w), 1 set
TH154-BW200	Wave jaws 40x200 mm (h x w), 1 set
	Inset jaw faces Quick change of jaw surface Price includes: 1 Set = 4 jaw faces
TH154-GT	Carrier for inset jaws, 1 set (4 carriers)
TH154-EB	Blank (smooth) inset jaws, 50x60 mm (h x w), 1 set
TH154-EBG	Rubber inset jaws, 50x60 mm (h x w), 1 set
TH154-EBP	Pyramid inset jaws, 50x60 mm (h x w), 1 set
TH154-EBV4	V-Inset jaws for round samples Ø 4 - 24 mm, 1 set
TH154-EBW	Wave inset jaws, 50x60 mm (h x w), 1 set
TH154-EBD	Diamond coated inset jaws, 50x60 mm (h x w), 1 set
TH154-GT100	Carrier for inset jaws, 1 set (4 carriers)
TH154-EB100	Blank (smooth) inset jaws, 50x100 mm (h x w), 1 set
TH154-EBG100	Rubber inset jaws, 50x100 mm (h x w), 1 set
TH154-EBP100	Pyramid inset jaws, 50x100 mm (h x w), 1 set
TH154-EBW100	Wave inset jaws, 50x100 mm (h x w), 1 set
TH154-EBD100	Diamond coated inset jaws, 50x100 mm (h x w), 1 set



www.grip.de/pdf/THS90_Datasheet.pdf


THS90



THS90	Vise grip, 50 kN, symmetric type of construction Price includes: 1 pair of vise grips
	Jaw faces for THS90 Price includes: 1 Set = 4 jaw faces
THS90-B	Blank jaws (smooth) without coating, 60x66 mm (h x w), 1 set
THS90-BG	Rubber jaws, 60x66 mm (h x w), 1 set
THS90-BP	Pyramid jaws (serrated), 60x66 mm (h x w), 1 set
THS90-BV	V-jaws for round samples Ø 5-30 mm, 60x66 mm (h x w), 1 set
THS90-BW	Wave jaws, 60x66 mm (h x w), 1 set
THS90-BD	Diamond coated jaws, 60x66 mm (h x w), 1 set
THS90-B210	Blank jaws (smooth) without coating, 70x210 mm (h x w), 1 set
THS90-BP210	Pyramid jaws (serrated), 70x210 mm (h x w), 1 set
THS90-BW210	Wave jaws, 70x210 mm (h x w), 1 set
	Further jaw surfaces and sizes on request



AML Instruments Ltd
 Vice Grips Data Sheet
 Contact us for friendly profesional advice Tel +44 (0) 1522 789375

		THS355
	THS355-200	Vise grip with 3 screws, 50 kN Body: aluminium Width: 200 mm Price includes: 1 pair of vise grips
	THS355-BW200	Wave jaws for textiles, 200 mm, 1 set
		THS668
	THS668-80	Vise grip, 100 kN Clamping force: to keep 200 kN with attraction force 1500 Nm using pneumatic screwdriver Price includes: 1 pair of vise grips
	THS668-80-BP	Pyramid jaws, 100 x 80 mm (h x w), 1 set
	THS668-210	Heavy vise grip with 8 M36 screws (pneumatic screwdriver recommended) Price includes: 1 pair of vise grips
	THS668-210-BP	Pyramid jaws, 70x210 mm (h x w), 1 set