

Manual Sample Cutting PRESSES “MP” series

(operation by hand lever)

Precision Strip Punches to cut RECTANGULAR, SQUARE or CIRCULAR precision samples with high accuracy, +/- 0,1mm, in order to evaluate *Tensile, basis weight, SCT, CMT, RCT, CCT...* for Paper (Liner / Medium) ...

According to Standards: ISO 3035 - ISO 7263 - ISO 3037 - ISO 536 - ISO 1924-2 - TAPPI T822 - TAPPI T809 - TAPPI T494...

- Easy to operate
- Rapidity and precision
- Robust Construction

RECTANGULAR versions:

- 15 x 300 mm **MP-15 model**
- 12,7 x 152,4 mm **MP-10 model**
- 25 x 300 mm **MP-25 model**
- Accuracy: +/- 0,1 mm



CIRCULAR versions:

- 100 cm² **MP-100 model**
- 50 cm² **MP-50 model**
- Accuracy: +/- 0,1 mm



Robust Strip Punches for preparation of test strips towards use in tensile strength, basis weight, Cobb, SCT, CMT, RCT, CCT, CLT and other similar tests.

Punch & Die manufactured from tempered tool alloy steel to ensure long life.

* TECHLAB SYSTEMS reserves the rights to make any technical modification without advance notice

Doc. : PRESSES_MP-1-CAT-I-R3