



Die Maschine

Rockwell & Brinell Hardness tester

Link
RB

Pro Series



The Link RB model is a Rockwell hardness tester that is capable of applying loads of 60, 100 and 150 Kg. It also has additional loads for Brinell tests (31.25, 62.5 y 187.5 Kg).

Loads are applied via a motor. Loads are easy to select using a rotary controller, and application speed is controlled via a hydraulic buffer.



Characteristics:

- Selectable Rockwell and Brinell loads.
- Load selection via rotary controller.
- Loads applied via motor.

Characteristics

| | |
|-------------------------------|----------------------|
| Preload Kg | 10 |
| Rockwell scale loads Kg | 60 - 100 - 150 |
| Brinell loads Kg | 31.25 - 62.5 - 187.5 |
| Vertical capacity mm | 170 |
| Goose neck (Horizontal) mm | 165 |
| Dimensions (W x D x H) mm | 215 x 520 x 700 |
| Net weight Kg | 100 |
| Power supply v | 220 |



Included accessories:

- Rockwell indenter - Ball 1/16".
- Rockwell indenter - Diamond cone.
- Brinell indenter - Ball Ø2.5 mm.
- Spare balls.
- Rockwell template (HRC, HRB).
- Horizontal support table Ø72 mm.
- Flat support table Ø150 mm (Optional).
- V-shaped support table Ø40 mm.
- Case with accessories.
- Instruction manual.
- Brinell conversion tables.



3 Year Warranty



with all Hoytom[®] machines

When you purchase a Hoytom[®] machine the equipment comes with a 3 year warranty.

Die Maschine



www.hoytom.com

 hoytom@hoytom.com



[youtube.com/hoytom](https://www.youtube.com/hoytom)



[hoytom](#)



[hoytom](#)