



# Impact and Thermal Printers



Model CBM-910



Model STP-103

## ***Model CBM-910 Palm-Sized Impact Printer***

### ***The Compact, Palm-Sized Printer with Unparalleled Features!***

The Model CBM-910 desk-top or wall mountable impact printer is light weight, ultra compact and perfect for a wide array of applications.

## ***Model STP-103 Thermal Printer***

### ***Speed, Quality, Quiet, Cost Effective***

The Model STP-103 is the ideal small footprint thermal printer. It is user-friendly and offers unsurpassed reliability.

#### ***Unique Features:***

- Compact and lightweight
- Easy paper loading
- Paper near end detection
- Included AC power adapter

# Impact and Thermal Printers

## Specifications

Model	CBM-910
Printing Method	Dot matrix
Power	DC 7V $\pm$ 1V (Printing), use exclusive adapter (DC 7V 1.6A)
Operating Temperature	30 - 104° F / 0 - 40° C
Dimensions (W x D x H)	4.17 x 7. x 3.5 in / 106 x 180 x 88 mm
Weight	1lb / 470 g
Safety Standards	CE Marking / TUV-GS / FCC Class, B / VCCI Class B, UL, C-UL (AC adapter PSE / U4 / UL / C-UL / KTL / CE / CCC / TUV-GS)

## Accessories

Paper..... 80251932

Printer Cable Kit..... 80252581

Pioneer™ Ranger™, Ranger™ Count

Cable Adapter\*\* ..... 80252583

Adapter..... 80251930

\*Printer Cable Kits include RS232 cable and cable adapter

\*\*Cable adapter connects Scout® Pro and Traveler™ models to the printer when the RS232 option module is installed

## Specifications

Model	STP-103
Printing Method	Direct thermal
Paper Sensor	Paper near end sensor, paper end sensor
Paper Type	Thermal paper
Power Supply	AC120 ~ 230V $\pm$ 10% (included)
Dimensions (W x L x H)	3.9 x 7.5 x 3.5 in / 100 x 191 x 90 mm
Safety Standards	CE Marking / TUV-GS / FCC Class B / VCCI Class B (AC adapter PSE / UL /C-UL / TUV-GS)

## Accessories

Thermal Paper. .... 80251931

Ink Ribbon..... 80251933

Printer Cable Kit 2\*. .... 80252571

Pioneer™, Ranger™, Ranger™ Count

Adapter..... 80251934

\*RS232 cable and the pin adapter are included in Printer cable kit

### OHAUS CORPORATION

Tel: 800.672.7722

973.377.9000

Fax: 973.944.7177

[www.ohaus.com](http://www.ohaus.com)

*The management system  
governing the manufacture  
of this product is  
ISO 9001:2015 certified.*

